

Nuclear Product Detection From Deuterated Nanoparticles Under Phonon Stimulation

Dr. John Dodaro

Aquarius Energy, CEO

Stanford University, Visiting Scholar

Washington, DC

September 8, 2023



U.S. DEPARTMENT OF
ENERGY

Public Release Slide





Project Title:

Nuclear Product Detection From Deuterated Nanoparticles Under Phonon Stimulation

PI: Matteo Cargnello (mcargn@stanford.edu)

Project Outcomes:

Demonstrate real-time production of helium and MeV charged particles from deuterated metallic nanoparticle alloys under resonant stimulation



Hypothesis

► Statement of hypothesis




- deuterated metallic nanoparticles become LENR-active by (1) creating active sites with small D-D separation, (2) exciting coherent vibrations, and (3) facilitating an energy transfer pathway that produces helium-4 and MeV-scale charged particles at reaction rates significantly amplified above of a bare D₂ molecule.

► Summary of pre, active, and post measurements

- Pre/Post TEM analysis of catalyst
- Independent Variables:
 - diode irradiation frequency: 800-1000nm (targeting 2–30 THz beat mixing)
 - substrate temperature: 0-500C (RTDs)
 - reactor pressure: 1-20 bar, with pressure transducers
 - catalyst composition: particle shape/dimensions & support (TEM), alloy/dopant ratios (EDX)
- Dependent Variables:
 - helium concentration: sub-100 ppb (high-res mass spec)
 - charged particle counts: micro-Watt flux (particle detector)
 - catalyst temperature: 0.01C (RTDs)


Experimental Setup – Reactor Cell

- ▶ Stainless steel 304 reactor; KBr window
- ▶ Macor block holds ceramic heater & substrate
- ▶ Cell volume 300cc; reduced for faster mass spec signal response
- ▶ Deuterium (99.999% purity)
 - isotopic enrichment 99.8
 - $H_2 < 100$, $HD < 3000$ ppmv




Experimental Setup – Catalyst substrate

- ▶ Alumina substrate holds catalyst layer
- ▶ SiC charged particle detector to catalyst ~5 micron photoresist



Experimental Setup – Solid state diodes


- ▶ IR diodes 800-1000nm (RPMC) with 1-6 Watt output
- ▶ Separated by 4"–8" from catalyst surface
- ▶ Frequency-tuned by Arroyo TECSOURCE power supply




Data Acquisition

Measurement	Recording Method	Settings	Latency	Storage Media
Helium, Nitrogen, Neon/Argon	Hidden DLS-20 quadrupole mass spec	events, counts/sec (bkgd 15 c/s)	10-30 sec	Hidden MasSoft
Temperature (catalyst, substrate)	RTDs (resistance)	0–500C	< 1 ms	National Instruments DAQ & Labview
Pressure	Pressure transducer (voltage)	1–20 bar	< 1 ms	National Instruments DAQ & Labview
Charged Particle Counts & Energy	Ortec pre-amp & radiation detector bias supply (current, voltage)	events, counts/sec	< 1 ms	National Instruments DAQ & Labview





Data Acquisition – Mass spectrometer (no catalyst)



Data Acquisition – Mass spectrometer (helium calibration)





Data Acquisition – Si detector counts & energy for calibration (Am241)



Modeling

- Alpha particle penetration calculation through SiC depletion layer
- Possible DFT & band structure calculations for phonon spectra




Initial Test Plan

- ▶ **Synthesis:** grow colloidal Pd-nanomaterials on stable supports
 - Carbon nanotube and ceramic supports for high thermal stability
 - Alloy particles to enhance non-linear optical response
- ▶ **Laser Stimulation:** Verify diode output at target frequency & THz mixing response
- ▶ **Particle Detector:** Work with UMich, TTU (Cat. B) teams for detector fabrication

Initial Results


► Developed Initial Catalyst Library

- Nanoparticles on CNT supports
- Alloy combinations, Pd-Ni
- Size 2-20nm, core-shell



► Demonstrated helium sensitivity

- Calibration pulse with air injection
- Two modes of operation: real-time “flow” runs & sustained “batch” runs



Plans for Next Quarter (Oct-Dec)

- ▶ Synthesize target nanoparticle material with reproducible recipe
 - Characterize with TEM, elemental analysis, chemisorption
 - Demonstrate target alloy/dopant concentrations
- ▶ Demonstrate tunable beat frequencies across 2–30 THz range
 - Explore nano-engineered approaches to maximizing laser mixing on catalyst (e.g. nanoparticle composition, dimensions, shape)
- ▶ Begin evaluation of SiC charged particle detector
 - Collect I-V curves demonstrating Schottky diode response
 - Record alpha counts & energy with source in reactor cell
- ▶ Demonstrate helium calibration sensitivity
 - Explore improvements, e.g. getter materials
 - Work with TTU team to test secondary validation